


MADE IN GERMANY · 47° 57' 31.4" N · 8° 46' 58.6" E



**GURBISCH**

## **Blepharo- plasty**

Instrumente für die Blepharoplastik  
Instruments for Blepharoplasty  
Instrumentos para la blefaroplastia  
Instruments pour la blépharoplastie  
Strumenti per la blefaroplastica


**medicon**<sup>®</sup>




**Instrumente für die Blepharoplastik**  
**Instruments for Blepharoplasty**  
**Instrumentos para la blefaroplastia**  
**Instrumentos pour la blépharoplastie**  
**Strumenti per la blefaroplastica**

Feine Wundhaken mit Isolierungsschicht  
 Fine retractors with insulation layer  
 Separadores finos con aislamiento  
 Ecarteurs fins avec revêtement isolant  
 Divaricatori fini con strato isolante


- VORTEILE**
  - ▶ Verhindert Kontakt zwischen HF-Nadel und Wundhaken
  - ▶ Vermeidet Gewebeschädigungen durch punktuelle Hitzeentwicklung
- ADVANTAGES**
  - ▶ Prevents contact between HF-needle and retractor
  - ▶ Avoids tissue damages caused by punctual heat development
- VENTAJAS**
  - ▶ Evita el contacto entre la aguja de alta frecuencia y el separador
  - ▶ Impide del traumatismo de tejido por el desarrollo de calor puntual
- AVANTAGES**
  - ▶ Empêche tout contact entre l'aiguille HF et l'écarteur
  - ▶ Evite le développement ponctuel de chaleur n'endommage les tissus
- VANTAGGI**
  - ▶ Impedisce il contatto tra l'ago ad alta frequenza e il divaricatore
  - ▶ Evita danni ai tessuti causati dal calore locale




**Fomon**



**Desmarres**



**Strandell-Stille**




**Langenbeck-Green**




**Instrumente für die Blepharoplastik**  
**Instruments for Blepharoplasty**  
**Instrumentos para la blefaroplastia**  
**Instrumentes pour la blépharoplastie**  
**Strumenti per la blefaroplastica**

Messzirkel  
 Caliper  
 Compás  
 Compas  
 Compasso



Blefaropinzette  
 Blepharoplasty forceps  
 Pinza para blefaroplastia  
 Pince de blépharoplastie  
 Pinza per blefaroplastica

Spezielle Biegung ermöglicht punktuellen Fassen des Gewebes  
 Special bending permits punctual tissue grabbing  
 Gracias al doblamiento especial el agarre de tejido puntual es posible  
 Sa courbure spéciale permet de saisir le tissu de manière ponctuelle  
 La curvatura speciale rende possibile la presa precisa del tessuto




Siebvorschlag  
 Set suggestion  
 Proposición de juego  
 Proposition pour boîte  
 Proposte di set

Item number	Description	Quantity
<b>04.21.14</b>	BROPHY (Sullivan) Delicate dissecting scissors · curved · sharp · 14 cm	1
<b>11.74.25</b>	BARRAQUER Micro needle holder · round handle Ø 8 mm · plain · 14 cm	1
<b>20.11.01</b>	FOMON Alar retractor · 17 cm	1
<b>20.11.02</b>	DESMARRES Lid retractor · 14 cm	1
<b>20.11.03</b>	STRANDELL STILLE Retractor · 6 x 18 mm · serrated · 19 cm	2
<b>20.11.04</b>	LANGENBECK-GREEN Retractor · 6 x 25 mm · serrated · 17 cm	1
<b>60.00.40</b>	CASTROVIEJO Caliper · 20 mm · 8 cm	1
<b>60.05.95</b>	JAEGER Lid plate · 11 cm	1
<b>60.33.03</b>	ANEL Tear duct cannulas · 23 Gauge · 4 cm · Luer-Lock	1
<b>60.33.17</b>	MOD. BONN Irrigation cannula · 19 Gauge · 3.5 cm	1
<b>60.34.00</b>	BOWMAN Lacrimal probe · figure 0000-000 · 13 cm	1
<b>60.34.21</b>	BOWMAN (Williams) Lacrimal probe · figure 00-0 · 13 cm	1
<b>60.34.81</b>	WILDER Lacrimal dilator · figure 1 · 11 cm	1
<b>60.34.82</b>	WILDER Lacrimal dilator · figure 2 · 11 cm	1
<b>60.47.05</b>	CASTROVIEJO Suturing forceps · 0,5 mm · 1 x 2 teeth · 10 cm	1
<b>60.55.45</b>	GUBISCH Blepharoplasty forceps · 10 cm	1
<b>60.73.93</b>	MOD. BONN Iris scissors · curved · 9 cm	1
<b>87.32.12</b>	Needle electrode · straight · Ø 4 mm	1
<b>87.33.12</b>	Needle electrode · straight · Ø 2.4 mm	1
<b>88.67.05</b>	Coagulation forceps bipolar · sharp · 16.5 cm	1
<b>49.50.10</b>	Medicontainer · 290 x 290 x 100 mm	1
<b>49.90.00</b>	Tray · 190 x 102 x 19 mm · with silicone mat	1
<b>49.97.90</b>	Wire basket · 240 x 250 x 40 mm · with silicone mat	1



Messbereich  
 Measuring range  
 Alcance de medición  
 Etendue de mesure  
 Campo di misura



Instrumente, optional  
 Instruments, optional  
 Instrumentos, opcional  
 Instruments, optionnel  
 Strumenti, opzionale

||| CMF

||| SURGICAL

||| NEURO+SPINE



**Medicon eG**

Gänsäcker 15  
D-78532 Tuttlingen

P. O. Box 44 55  
D-78509 Tuttlingen

Tel.: +49 (0) 74 62 / 20 09-0  
Fax: +49 (0) 74 62 / 20 09-50

E-Mail: [sales@medicon.de](mailto:sales@medicon.de)  
Internet: [www.medicon.de](http://www.medicon.de)

Germany

Modelländerungen vorbehalten  
Patterns are subject to change  
Salvo modificaciones  
Tous droits réservés des changements de modèle  
Ci riserviamo la facoltà di cambiamenti nei modelli

Gedruckt in Deutschland  
Printed in Germany  
Impresso en Alemania  
Imprimé en Allemagne  
Stampato in R.F.G.  
Trykkes i Tyskland

© Copyright 03/2013

451 .02.05

**medicon**<sup>®</sup>

